



the Dropper Plus

Point-of-Care Urine Dipstick Control

- Liquid, human urine-based
- **1 month open vial stability at room temperature for maximum portability**
- 18 month unopened and open vial stability when stored at 2-8 C*
- Quantrol™ Quality Control Program available

Analytes*

Bilirubin, Urine	Glucose	Leukocytes	pH	Total Protein
Blood	hCG	Microalbumin	Protein	Urobilinogen
Creatinine	Ketones	Nitrite	Specific Gravity	

**Refer to the package insert of currently available lots for specific analyte and stability claims*

The Dropper Plus is specifically designed for point-of-care dipstick control. Its unique, one month room temperature stability eliminates storage problems and provides maximum portability – ideal for sites without refrigeration. A new feature is the addition of creatinine and microalbumin. Values are provided for all analytes; qualitative results for hCG for early pregnancy detection test methods are included. This control is intended for use with Multistix[®], Chemstrip[®], DiaScreen[®], UriScan™, pH (Nitrazine paper) and other urinalysis reagent strips. It can also be used for confirmatory tests such as Acetest[®], Clinitest[®], Ictotest[®] reagent tablets and refractometry.

ORDERING INFORMATION:

1440-04 2 Level Set.....10 x 5 mL



For More Information: phone +1.310.536.0006 email info@4qc.com web www.4qc.com