

UrichemTRAK®

Liquid Stable Urine Chemistry Control

UrichemTRAK® is a liquid ready to use control designed to monitor methods for routine urine chemistry analytes.

Key Features

- Liquid, ready to use; no preparation required
- Human urine based
- Non-azide formulation will not interfere with ISEs
- Includes general chemistry, physical properties and hCG tests
- LabLink® xL QAP available for all users
- Product stability: *
 - 2 year opened **and** unopened-vial stability at 2-8°C



Analytes*

- | | | | |
|------------|--------------|------------|------------------|
| Amylase | Glucose | pH | Specific Gravity |
| BUN | hCG | Phosphorus | Total Protein |
| Calcium | Magnesium | Potassium | Uric Acid |
| Chloride | Microalbumin | Sodium | |
| Creatinine | Osmolality | | |

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
UrichemTRAK®	Level 1	6 x 15 mL	UR11001
	Level 2	6 x 15 mL	UR22002
	Multi-Pack,	6 x 15 mL	UR-MP
	Bi-Level		