

Liquid Stable Immunoassay Control

Liquimmune® is a comprehensive immunoassay control for routine immunoassay parameters on a full range of immunoassay platforms.

Key Features

- Liquid, ready to use; no preparation required
- 3 significant levels
- Human serum based
- LabLink® xL QAP available for all users
- Product stability:*
 - 30 day open-vial stability at 2-8°C for all analytes
 - 90 day unopened-vial stability at 2-8°C
 - 3 year shelf life at -20°C

Analytes*

17- α -hydroxyprogesterone	Estradiol	Lidocaine	T-Uptake
Acetaminophen	Ethosuximide	Lithium	T3
Aldosterone	Ferritin	Luteinizing Hormone (LH)	T3 (Free)
alpha-Fetoprotein (AFP)	Folate	NAPA	T4
Amikacin	FSH	Phenobarbital	T4 (Free)
Benzodiazepines	Gentamicin	Phenytoin	Testosterone
beta-hCG	hCG	Primidone	Theophylline
Carbamazepine	Homocysteine	Procainamide	Thyroglobulin
CEA	Human Growth Hormone	Progesterone	Tobramycin
Chloramphenicol	Immunoglobulin A	Prolactin	Tricyclic Antidepressants
CK-MB	Immunoglobulin E	PSA	TSH
Cortisol	Immunoglobulin G	PSA (Free)	Valproic Acid
DHEA-Sulfate	Immunoglobulin M	PTH**	Vancomycin
Digoxin	Insulin	Quinidine	Vitamin B12
Disopyramide	LD-1 Isoenzyme	Salicylate	

* Refer to package insert for specific analyte and stability claims
 ** Coming soon

Ordering Information

Product Description	Level	Pkg. Size	Product Code
Liquimmune®	Level 1	6 x 5 mL	LIG-101
	Level 2	6 x 5 mL	LIG-202
	Level 3	6 x 5 mL	LIG-303
Dropper Tips		Pkg. 100	286-606

