

Liquid Stable Unassayed Therapeutic Drug Monitoring Control

chemTRAK® TDM is an unassayed, multi-constituent control containing widely tested TDMs and immunoassay parameters.

Key Features

- Liquid, ready to use; no preparation required
- 3 significant levels
- LabLink® xL QAP available for all users
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at 2-8°C



Analytes*

- | | | | |
|----------------------|---------------|------------------|---------------------------|
| Acetaminophen | Ethanol | Phenytoin (Free) | Tobramycin |
| Amikacin | Ethosuximide | Primidone | Tricyclic Antidepressants |
| Amitriptyline | Gentamicin | Procainamide | TSH |
| Caffeine | Lidocaine | Quinidine | Valproic Acid |
| Carbamazepine | Lithium | Salicylate | Valproic Acid (Free) |
| Carbamazepine (Free) | Methotrexate | T3 | Vancomycin |
| Chloramphenicol | NAPA | T4 | |
| Cyclosporine | Netilmicin | Theophylline | |
| Digoxin | Nortriptyline | | |
| Disopyramide | Phenobarbital | | |
| Estriol | Phenytoin | | |

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
chemTRAK® TDM	Level 1	6 x 5 mL	TD1-1001
	Level 2	6 x 5 mL	TD2-2002
	Level 3	6 x 5 mL	TD3-3003
Dropper Tips		Pkg. 100	286-606

