

Liquid Stable Cardiac Marker Control

CardioImmune®•TL is a comprehensive Cardiac Marker Control for monitoring common and newer Cardiac Marker assays and is assayed for all major platforms.

Key Features

- Liquid, ready to use; no preparation required
- 4 distinct levels
- Ultra-low levels in Level L
- Human serum based
- LabLink® xL QAP available for all users
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at -20°C



Coming in Summer 2005 NEW
 CardioImmune®•xL with the addition of BNP-32 & Homocysteine

Analytes*

- | | | |
|-----------|-----------|------------|
| hsCRP | Myoglobin | Troponin I |
| CK-MB | ProBNP | Troponin T |
| Digitoxin | | |

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
CardioImmune®•TL	Level L	6 x 3 mL	CAI-TLL
	Level 1	6 x 3 mL	CAI-TL1
	Level 2	6 x 3 mL	CAI-TL2
	Level 3	6 x 3 mL	CAI-TL3
	Multi-Pack, Tri-Level	6 x 3 mL	CAI-TL4
Dropper Tips		Pkg. 100	286-606

