



Liquid Stable Diabetes Control

MAS Liquid Diabetes Control is a convenient ready to use whole blood control for use on automated HPLC and Immunoassay HbA_{1c} instruments.

Key Features

- Liquid, ready to use; no preparation required
- Human whole blood based
- Convenient multi-pack format
- Aliquot and re-freeze supported
- Assayed for most major instrument systems
- A true 3rd party assessment of Bio-Rad Variant[®] / Variant[®] II / D10[®] systems
- LabLink[®] xL QAP available for all users
- Product stability:*
 - 14 day open-vial stability at 2-8°C
 - 2 year shelf life at -20°C



Analytes*

Hemoglobin A_{1c}

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS [®] Liquid Diabetes Control	Multi-Pack, Bi-Level	6 x 1 mL	DBCL-MP