

MAS[®] DOA-XAB[®]

Liquid Stable Drugs of Abuse Control

MAS[®] DOA-XAB[®] is designed to monitor drugs of abuse screening test methods for the Abbott[®] reagent systems.

Key Features

- Liquid, ready to use; no preparation required
- Human urine based
- Contains the SAMHSA** 5, plus Barbiturates, Benzodiazepines, Ethanol, LSD, Methadone, Methaqualone and Propoxyphene
- Cut-offs are in accordance with SAMHSA and leading manufacturers' guidelines
- Available in four distinct levels
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at 2-8°C



Analytes*

Benzoylcegonine
Cannabinoid (11-nor-9-COOH- Δ^9 -THC)
d-Methamphetamine

Ethanol
LSD
Methadone
Methaqualone

Morphine
Nordiazepam

Phencyclidine
Propoxyphene
Secobarbital

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS[®] DOA-XAB[®]	Level 1, negative	6 x 18 mL	DOA-101
	Level 2, 20-25% below challenge to immunoassay cut-off	6 x 18 mL	DOA-202
	Level 3, 20-25% above challenge to immunoassay cut-off	6 x 18 mL	DOA-303
	Level 4, elevated immunoassay	6 x 18 mL	DOA-404

** Substance Abuse and Mental Health Services Administration (SAMHSA)

