

Liquid Stable Assayed and Unassayed Chemistry Control

chemTRAK® Platinum is a comprehensive tri-level liquid general chemistry and TDM control.

Key Features

- Liquid, ready to use; no preparation required
- A complete menu of general chemistry and TDM analytes
- 3 levels available (low; mid; high)
- Assayed values available for major platforms
- LabLink® xL QAP available for all users
- Product stability:*
 - 14 day open-vial stability at 2-8°C
 - 60 day unopened-vial stability at 2-8°C
 - 2 year shelf life at -20°C

Analytes*

Acetaminophen	Carbamazepine	Imipramine	Primidone
Acid Phosphatase (ACP)	Chloramphenicol	Iron	Procainamide
Albumin (BCG)	Chloride	Iron Binding Capacity	Quinidine
Alkaline Phosphatase (ALP)	Cholesterol	Lactic Acid	Salicylate
ALT (GPT)	Creatine Kinase (CK)	LDH	Sodium
Amikacin	CO ₂	Lidocaine	T-Uptake
Amylase	Creatinine	Lipase	T4
AST (GOT)	Cyclosporine	Lithium	Theophylline
Bilirubin, Conjugated	Digoxin	Magnesium	TIBC
Bilirubin, Direct	Disopyramide	Methotrexate	Tobramycin
Bilirubin, Neonatal	Estriol	NAPA	Total Protein
Bilirubin, Total	Ethanol	Netilmicin	Tricyclic Antidepressants
Bilirubin, Unconjugated	Ethosuximide	Phenobarbital	Triglycerides
BUN	Gentamicin	Phenytoin	Uric Acid
Caffeine	GGT	Phosphorus	Valproic Acid
Calcium	Glucose	Potassium	Vancorycin

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
chemTRAK® Platinum (Assayed)	Level 1	10 x 15 mL	AP11001
	Level 2	10 x 15 mL	AP22002
	Level 3	10 x 15 mL	AP33003
chemTRAK® Platinum (Unassayed)	Level 1	10 x 15 mL	UP11001
	Level 2	10 x 15 mL	UP22002
	Level 3	10 x 15 mL	UP33003
Dropper Tips		Pkg. 100	286-606

