

MAS[®] DOA GC/MS

Liquid Stable Drugs of Abuse Control

MAS[®] DOA GC/MS is designed to monitor the quantitative confirmatory test methods of gas chromatography and mass spectrometry systems.

Key Features

- Ready to use; no preparation required
- Human urine based
- Cut-offs are in accordance with SAMHSA** guidelines
- Available in four distinct levels
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at 2-8°C



Analytes*

Amobarbital	Codeine	Methadone	Pentobarbital
Benzoylcegonine	d-Amphetamine	Methaqualone	Phencyclidine
Butalbital	d-Methamphetamine	Morphine-3-β-d-	Phenobarbital
Cannabinoid (11-nor-9-COOH-Δ ⁹ -THC)	Ethanol	Glucuronide	Propoxyphene
	LSD	Nordiazepam	Secobarbital

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS [®] DOA GC/MS	Level 1, negative	6 x 18 mL	DOA-101
	Level G2, 20-25% below GC/MS cut-off	6 x 18 mL	DOA-505
	Level G3, 20-25% above GC/MS cut-off	6 x 18 mL	DOA-606
	Level G4, elevated GC/MS	6 x 18 mL	DOA-707