

MAS[®] Opiates 2000

Liquid Stable Opiates DOA Control

MAS[®] Opiates 2000 is designed to monitor semi-quantitative Opiates screening methods at the SAMHSA** cut-off of 2000 ng/mL.

Key Features

- Liquid, ready to use; no preparation required
- Human urine based
- Available in two levels; 25% below and 25% above the 2000 ng/mL cut-off
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at 2-8°C



Analytes*

Morphine

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS [®] Opiates 2000	Multi-Pack, Bi-Level	4 x 18 mL	OP2000

