

MAS[®] Opiates 2000 GC/MS

Liquid Stable Opiates DOA Control

MAS[®] Opiates 2000 GC/MS is designed to monitor quantitative and confirmatory Opiate test methods utilising Opiates screening methods at the SAMHSA** cut-off of 2000 ng/mL.

Key Features

- Liquid, ready to use; no preparation required
- Human urine based
- Available in two levels; 25% below and 25% above the 2000 ng/mL cut-off
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 3 year shelf life at 2-8°C



Analytes*

6-Acetylmorphine

Codeine

Morphine-3-β-d-Glucuronide

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description

MAS[®] Opiates 2000 GC/MS

Level

Multi-Pack, Bi-Level

Pkg. Size

4 x 18 mL

Product Code

OP2GCMS