


MAS[®] Alcohol/Ammonia

Liquid Stable Assayed Alcohol/Ammonia Control

MAS[®] Alcohol/Ammonia is a liquid, assayed, bi-level control for monitoring alcohol and ammonia levels in serum.

Key Features

- Liquid, ready to use; no preparation required
- Human serum based
- Normal ammonia values in Level 1
- Convenient multi-pack format
- Assayed values available for major instruments
- LabLink[®] xL QAP available for all users
- Product stability:*
 - 30 day open-vial stability at 2-8°C
 - 2 year shelf life at 2-8°C



Analytes*

Ethanol

Ammonia

* Refer to package insert for specific analyte and stability claims

Ordering Information

| Product Description | Level | Pkg. Size | Product Code |
|----------------------------------|-------------------------|------------|--------------|
| MAS [®] Alcohol/Ammonia | Multi-Pack, Bi-Level | 6 x 3.5 mL | AAC-MP |
| Dropper Tips | | Pkg. 100 | 286-606 |

