

## MAS<sup>®</sup> Bilirubin

### Liquid Stable Assayed Bilirubin Control

MAS<sup>®</sup> Bilirubin is an assayed, tri-level control designed to monitor the performance of total, conjugated, unconjugated and neonatal bilirubin.

#### Key Features

- Liquid, ready to use; no preparation required
- Extremely stable, complements any General Chemistry Control
- 3 levels offered to assess the performance of the full assay range
- Elevated bilirubin in Level 3 for neonatal applications
- LabLink<sup>®</sup> xL QAP available for all users
- Product stability:\*
  - 14 day open-vial stability when stored in the dark at 2-8°C
  - 2 year shelf life when stored in the dark at 2-8°C



#### Analytes\*

Bilirubin, Conjugated  
Bilirubin, Direct

Bilirubin, Neonatal  
Bilirubin, Total

Bilirubin, Unconjugated

\* Refer to package insert for specific analyte and stability claims

#### Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS <sup>®</sup> Bilirubin	Level 1	6 x 5 mL	BC101
	Level 2	6 x 5 mL	BC102
	Level 3	6 x 5 mL	BC103
Dropper Tips		Pkg. 100	286-606