

Liquid Stable Assayed Cerebrospinal Fluid Control

MAS[®] CSF is an assayed, bi-level control designed to monitor the performance of quantitative and manual chemical analysis of human cerebrospinal fluid.

Key Features

- Liquid, ready to use; no preparation required
- Oligoclonal bands available in Level 2
- Assayed for cellulose acetate and agarose gel electrophoresis
- LabLink[®] xL QAP available for all users
- Convenient multi-pack format
- Product stability:*
 - 30 days open-vial stability at 2-8°C
 - 2 year shelf life at 2-8°C



Analytes*

Chloride	Immunoglobulin G	Oligoclonal Bands	Total Protein
Globulins	Immunoglobulin M	Potassium	
Glucose	Lactic Acid	Prealbumin	
Glutamine	LDH	Protein Electrophoresis	
Immunoglobulin A	Microalbumin	Sodium	

* Refer to package insert for specific analyte and stability claims

Ordering Information

Product Description	Level	Pkg. Size	Product Code
MAS [®] CSF	Multi-Pack, Bi-Level	6 x 3.5 mL	CSF-MP
Dropper Tips		Pkg. 100	286-606

