

**Liquid Stable Lipid Control**

MAS<sup>®</sup> Lipid is an assayed, bi-level control offering a complete lipid profile.

Serum Chemistry Controls

**Key Features**

- Liquid, ready to use; no preparation required
- Human serum based
- No product freezing or aliquoting required
- Analyte levels are at critical decision points for your total lipid profile
- Assayed for direct HDL and LDL methodologies
- LabLink<sup>®</sup> xL QAP available for all users
- High Sensitivity CRP and Lipoprotein (a) included
- Product stability:\*
  - 30 days open-vial stability at 2-8°C
  - 2 year shelf life at 2-8°C



**Analytes\***

Apolipoprotein A  
Apolipoprotein B  
Cholesterol

HDL Cholesterol  
hsCRP  
LDL Cholesterol

Lipoprotein (a)  
Phospholipids  
Triglycerides

\* Refer to package insert for specific analyte and stability claims

**Ordering Information**

**Product Description**

**MAS<sup>®</sup> Lipid**

**Dropper Tips**

**Level**

Level 1  
Level 2

**Pkg. Size**

6 x 3 mL  
6 x 3 mL  
Pkg. 100

**Product Code**

LIP-101  
LIP-202  
286-606