

VITREX



Comprehensive range of high quality capillary blood lancets and lancing devices

- ◆ LOW COST LANCETS - for general clinical use
- ◆ LANCETS - for blood glucose monitoring
- ◆ SINGLE USE LANCETS - for all multi-patient and clinical laboratory use
- ◆ VARIABLE DEPTH LANCING DEVICES

VITREX

Distributed by:
Vitrex Medical Limited
Manorfields House, Duncote,
Towcester, Northants NN12 8AH
United Kingdom
T +44 1327 358157
F +44 1327 350421
E vitrex@vitmed.com

Vitrex Medical Aps
P.O.B. 507
Vasekaer 6-8
2730 Herlev
Denmark

CE 0459

01/05/03
Edition 2



VITREX STEEL

▶ GENERAL PURPOSE LANCETS

Vitrex Steel: Hand held, flat, stainless steel blood lancets
 Hardened surgical steel - rigidity and safety
 High quality, asymmetric point - excellent blood flow
 Sterile - gamma irradiation to DIN EN 552/556
 Convenient - 200 piece box
 Article Number: 3502



VITREX SOFT/GENTLE

▶ BLOOD GLUCOSE LANCETS

Vitrex Soft: 23g needle with plastic body
 Precision ground, 3 facet point - patient comfort with good blood flow
 Article Number: 3512

Vitrex Gentle: 28g needle with plastic body
 Silicone treated needle - reduced penetration force
 High quality, ultra sharp point - improved patient comfort
 Suitable for use with modern low sample volume meters
 Article Number: 3532

Common features
 Standard size - fit most commonly used lancing devices (not SOFTGLIX*)
 Sterile - gamma irradiation at min 25 kGy
 Convenient - 200 piece box



VITREX PREMIER/COMPACT

▶ LANCING DEVICES

Vitrex Premier: Standard lancing device
 Non-touch lancet removal - protects from accidental needlesick injury
 9 position adjustable depth setting - for all types and thicknesses of skin
 Article Number: 3690

Vitrex Compact: Compact lancing device
 5 position adjustable depth setting - for a variety of skin thicknesses
 Article Number: 3691

Common features
 Used to take samples for blood glucose testing
 Accept Vitrex and most other popular lancet brands
 Simple convenient activation and reactivation
 Individually packaged



VITREX SAFE

▶ SAFETY LANCETS

Vitrex Safe - Blue
 Penetration depth 1.8mm - for general use.
 Ultra sharp, high quality needle point - rapid, pain-free entry
 Article Number: 3582

Vitrex Safe - Green
 Penetration depth 2.4mm - for thicker skin
 Ultra sharp, high quality, needle point - rapid pain-free entry
 Article Number: 3592

Vitrex Safe - Yellow
 Penetration depth 1.0mm - neo-natal heelsticks/paediatric use
 Ultra sharp, high quality blade - increased blood flow
 Article Number: 3552

Common features
 Pressure activated release - consistently good blood samples
 Fully shielded needle - protects PHLEBOTOMISTS from injury
 Guaranteed single use - protects PATIENTS from cross-infection
 Silicone treated needle/blade - improves patient comfort
 Sterile - gamma irradiation at min 25 kGy
 Convenient - 200 piece boxes